


Price List

*Precision, Quality, and Value
Your Trusted Partner in Dental Restoration*

20
26



 info@mag-dental.com

 1300 538 992

 www.mag-dental.com



Crown & Bridge

All Ceramic

Zirconia – Monolithic	\$190
Zirconia – Layered	\$220
e.max® – Monolithic	\$220
e.max® – Inlay/Onlay	\$220
e.max® – Veneer	\$220
e.max® – Layered	\$280



PFM

PFM – Non-Precious	\$180
Full-Cast Metal Crown	\$210
180° Porcelain Margin	\$20
360° Porcelain Margin	\$35

Precious alloys additional at market price



Add Ons

Temporary PMMA Crown	\$80
Pink Porcelain (per unit)	\$40
Diagnostic Wax Up (per tooth)	\$35
3D Printed Full Models (per arch)	\$20
Putty Key for Temporary	\$35



Implant Solutions

Implant Package

All-Inclusive Package \$649
(package contains: Analog, screw, cement or screw retained crown (Zirconia/PFM), custom Ti-abutment and soft tissue model)



Full Arch Implant (per arch)

Full Arch Zirconia Bridge with \$3990
Titanium Frame (All-On-X)

Add Ons

Simple Implant Positioning Guide \$100
(includes one implant)

Additional Implant Position on Simple Surgical Guide \$45

Temporary Implant Crown \$100





Removable Appliances



Orthodontics

Bleaching Tray	\$80
Clear Retainer	\$90
Hawley Retainer	\$120
Space Maintainer	\$140
Twin Block	\$330
Bionator	\$330
Lingual Side Wire	\$90
Expander Appliances	\$300
Frankel Appliances	\$390
Custom Orthodontic Appliance	POR



Mouth Guards

Nylon Splint	\$295
Hard Night Guard	\$140
Soft Night Guard	\$140
Dual-Layer Night Guard (Soft/Hard)	\$160
Sports Mouth Guard (Single-Layer)	\$99
Sports Mouth Guard (Double-Layer)	\$110
Professional Sports Mouth Guard	\$150
Add Name on Mouth Guard	\$12



Snoring Devices

Connector Bar Anti-Snoring Appliance	\$400
Shark Fin Adjustable Apnea Appliance	\$500

Orthodontics Appliances

Expansion Appliances



Fixed Bidirectional
Dental Palate Expander



Removable
Bidirectional Teeth
Palate Expander



SS Nance Appliance




Stainless Steel Nance
Appliance Palatal Bar


Functional Appliances



Twin Block



Pendulum Appliance




Frankel III Malocclusion
Orthodontic Functional
Appliance



Synthetic Activator
Appliance

Retentive Appliances




Transparent Erkodont Silensor



Clear Retainer



Hawley Retainer



Space Maintainer



Lingual Wire Retaine

Scanner Setup Guide

1. 3Shape TRIOS (via 3Shape Communicate)

1. Log into your 3Shape Communicate account (communicate.3shape.com)
2. Go to the "Connections" tab
3. Click "Add Connection"
4. Enter our lab email: **info@mag-dental.com**
5. Click "Connect" and wait for confirmation
6. Once accepted, you can send scans to us directly

2. iTero Scanner (Align Technology)

1. Log into your myiTero account
2. Click the Settings icon → Select Practice Settings
3. Go to the My Labs tab → Click Add Lab
4. Search for "**MAG Dental Laboratory**" and click Connect
5. Confirm the connection — it will appear under My Labs once linked
6. After connection, select MAG Dental Laboratory when submitting new cases

3. Medit Scanner (via Medit Link)

1. Log into Medit Link
2. Go to the "Partners" tab
3. Click "Add Partner"
4. Search for: MAG Dental Laboratory or use email: **info@mag-dental.com**
5. Click "Request" to connect. Once approved, you can send scans directly

4. Primescan (via DS Core Platform)

1. Log in to your DS Core account
2. In the top-right corner, click on your account name and select "Practice"
3. From the Practice menu, go to "Preferred Labs"
4. Click "Add Preferred Labs"
5. In the search bar, type: **MAG Dental Laboratory**
6. The lab list will adjust as you type. Once you see our name, click "Add Lab"
7. A confirmation window will appear — click "Add Lab" again to confirm
8. Use the left arrow in the top-left corner to return to your Preferred Labs list
9. MAG Dental Laboratory will now appear with the status "Active"

Need Help?

If you're unsure how to complete any of the steps above, or your scanner brand is not listed, just reach out — we'll guide you through it or register with your system if needed.